

Intended Use:

RF (Rheumatoid Factor) Test reagent/kit is a medical device intended for the screening of Rheumatoid Factor (RF) in human specimen. It is intended for *in vitro* diagnostics use by trained healthcare professionals in IVD laboratory only and not near patients.

Clinical Significance

Measurement of rheumatoid factor is used for differentiating rheumatoid arthritis from other chronic inflammatory arthritis and is important in the progress and therapeutic management of the disease. Rheumatoid factor has been associated with some bacterial and viral infections (ex. Hepatitis, infectious, Mononucleosis) some chronic infections (ex. Tuberculosis, parasitic disease. Subacute Bacterial Endocarditis) and cancer.

Principle

The RF-Turbi latex is a quantitative turbidimetric test for the measurement of RF in human serum or plasma. Latex particles coated with human gamma globulin are agglutinated when mixed with samples containing RF. The agglutination causes an absorbance change, dependent upon the RF contents of sample that can be quantified by comparison from a calibrator of known RF concentration.

Reagent Composition

R1 Diluent	Tris buffer 20 mmol/L, pH 8.2.
R2 Latex	Latex particles coated with human gamma globulin, pH 8.2
RF Calibrator	Human serum RF concentration is stated on vial label. and it is traceable to the "Rheumatoid Arthritis Serum" 64/002 from WHO (NIBSC)

Precaution

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

Calibration

Use RF Calibrator ONLY.

The sensitivity of the assay and the target value of the calibrator have been

Standardized against the International Reference Standard from NIBSC 64/002. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

Preparation

RF Calibrator: Ready to use.

Calibration Curve: Prepare the following RF calibrator dilutions in NaCl 9g/L. Multiply the concentration of the RF calibrator by the corresponding factor stated in table below to obtain the RF concentration for each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator RF (μL)	--	10	20	40	60	80
NaCl 9g/L (μL)	80	70	60	40	20	--
Factor	0	0.125	0.25	0.5	0.75	1.0

Working Reagent:

Swirl the latex vial before use. Mix Latex and Diluent in a 1:4 ratio (i.e 2 ml R2 + 8 ml R1) prior to use.

Working Reagent is stable during 20 days at 2-8°C

Stability and Storage

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: Presence of particles and turbidity.

Additional Equipment

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 540 nm filter.

Samples

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

Assay Parameters

Mode	Two Point / Fixed Time
Reaction	Ascending
Wavelength	650 nm (600-650)
Blank With	Distilled Water
Sample Volume	5 μl
Reagent Volume	400 μl : 100 μl
Delay Time	10 (Sec)
Read Time	120 (Sec)
Calibrator Conc.	See Vial Label
Linearity Limit	Up to 100 IU/ml
Unit	IU/ml

Procedure

1. Bring the reagents and the photometer (cuvette holder) to 37°C.
2. Assay conditions:
Wavelength: 650 nm (600-650)
Temperature: 37°C
Cuvette light path: 1 cm
3. Adjust the instrument to zero with distilled water.
4. Pipette into a cuvette

Addition Sequence	Test
R1 Diluent	400 μl
R2 Latex	100 μl
Calibrator or Sample	5 μl

5. Mix and read the absorbance immediately Abs(A1) and after 2 minutes Abs(A2) of the sample addition.

Calculation

$$\text{RF (IU/ml)} = \frac{\text{Abs (A2-A1) sample}}{\text{Abs (A2-A1) calibrator}} \times \text{Calibrator conc.}$$



OR

Calculate the absorbance difference

(A2-A1 blank reagent) of each point of the calibration curve and plot the values obtained against the RF conc. Of each calibrator dilution.

Rheumatoid Factor concentration in the sample is calculated by interpolation of its (A2-A1 blank reagent) in the calibration curve.

Quality Control

Control Sera are recommended to monitor the performance of manual and automated assay procedures.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

Reference Values

Normal values up to 30 IU/mL.

Each laboratory should establish its own reference range.

Performance Characteristics

1. Limit detection: Values less than 5 IU/ml give non-reproducible results.
2. Measurement range: 6-100 IU/ml, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9g/L and retested again.
3. The linearity limit and measurement range depends on the sample reagent/ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
4. Prozone effect: No prozone effect was detected up to 800 IU/ml.
5. Sensitivity: 3.0 mA.U./ml.
6. Precision: The reagent has been tested for 20 days, using three different RF concentrations in a study.

Notes

1. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

References

1. Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951- 960.
2. Robert W Dorner et al. Clinica Chimica Acta 1987; 167: 1-21.
3. Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528- 534.
4. Vladimir Muié et al. Scand J Rheumatology 1972; 1: 181- 187.
5. Paul R et al. Clin Chem 1979; 25/11: 1909- 1914.
6. Young DS. Effects of drug on clinical laboratory test, 4th ed. AACC Press, 1995

Packaging

R1. Diluent: 1 x 40 ml

R2. Latex : 1 x 10 ml

RF-CAL: 1 x 2 ml

For *in vitro* Diagnostic Use Only